Directions for Use and Infection Control Procedures

sylc® CeramicTip

Sylc® CeramicTip E-Type Air Polishing Device

DIRECTIONS FOR USE
The Sylc® CeramicTip E-Type Air Polishing Tip is compatible for use with EMS® AIR-FLOW® Handy 2+ air polishing handpiece for the removal of extrinsic stain, dental plaque and soft debris while cleaning, and polishing the tooth surface.

INDICATIONS: The Sylc® CeramicTip E-Type is specifically designed to dispense Sylc® Therapeutic Prophyl powder. The device may also be used with sodium bicarbonate and other air polishing powders.

DESCRIPTION: The Sylc® CeramicTip E-Type is constructed from wear resistant ceramic materials and medical grade steel.

CONTRAINDICATIONS: Do not direct the air polishing stream at soft tissue or into the sulcus.

PRECAUTIONS:

- During the preparation, process and application, it is recommended that both clinician and patient wear protective eyewear.
- If any material gets into the eyes, immediately flush with water. Do not direct the jet of powder toward fillings, or crown and bridgework as this could damage the restorations.
- A dampened gauze pad can be placed on the patient's lip to protect the skin and increase patient comfort.
- Delivering Sylc or other polishing powders with air polishing equipment at high air pressures (bar/psi), may cause injury to the gum tissues or an emphysema caused by the introduction of air into the soft tissue spaces. It is strongly recommended not to direct Sylc at the gum tissue or into the gingival sulcus.
- When using Sylc powder, please refer to the Sylc Therapeutic Prophyl Powder DFU for specific precautions.
- End of life of the Sylc CeramicTip E-Type is normally determined by routine use and wear.
- The Sylc CeramicTip E-Type lifetime will be optimized when working no higher than the recommended operating air pressure range of 2.8-3.2 bar (40-46 psi).
- Do not drop the Sylc CeramicTip E-Type as it may cause damage to the tip.
- Denfotex Research and the distributor of this device accept no fault for direct or consequential injury or damage resulting from improper use, arising in particular through non-observance of the operating instructions, or improper preparation and maintenance.

INSTALLATION: Follow the equipment manufacturer's guidelines for installing the Sylc CeramicTip E-Type to the EMS® AIR-FLOW® Handy 2+ handpiece. **Ensure that all condensed water is removed from the compressor and the air line is dry.** Slide the CeramicTip E-Type into the handpiece and gently push until fully seated. To remove, rotate the CeramicTip E-Type whilst pulling at a slight angle. Do not pull straight out.

USAGE:

- Before each use, disinfect, clean and sterilize the different pieces and accessories of this device.
- Ensure Sylc CeramicTip E-Type is totally dry prior to use to prevent powder build-up.
- If moisture is present in the Sylc CeramicTip E-Type, blow the interior of the tip with compressed air to avoid the formation of powder concretions in the air tubes.
- Install the Sylc CeramicTip E-Type onto EMS handpiece and run dry air through device before filling powder chamber.
- Apply the Sylc powder according to the Sylc Therapeutic Prophyl Powder directions for use Instructions.

sylc® CeramicTip

MANUFACTURER:
Denfotex Research Ltd.
Suite 2, Forth House
Burnside Business Court,
North Road, Inverkeithing
KY11 1NZ
United Kingdom

DENFOTEX

www.denfotexresearch.com

Patent Pending.
Sylc® is a registered trademark of Denfotex Research Ltd.
Electro Medical Systems S.A. has not authorized the Sylc® CeramicTip E-Type design.
Air-Flow® is a registered trademark of EMS- Electro Medical Systems S.A.

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Sylc® è un marchio registrato di Denfotex Research Ltd.
Electro Medical Systems S.A. non ha autorizzato il modello Sylc® CeramicTip E-Type.
Air-Flow® è un marchio registrato di EMS-Electro Medical Systems S.A.